



ICOS



INTEGRATED
CARBON
OBSERVATION
SYSTEM



EcoTech FTIR instrument modifications and testing for tall tower deployment

Alex Vermeulen, Arjan Hensen, Pim van den Bulk,
Chris Edwards, Marie Laborde, David Griffith



Instrumentation needs to be:

- Robust
- Very high precision
- High time resolution
- Continuous measurement
- Lifetime >10 years
- Low maintenance
- Low consumables
- Easy to operate
- Completely automated
- Easy to integrate
- Low power+gas consumption
- Commercially available
- Supplier support



Tall tower applications

- Vertical gradients along tower
 - Switching levels
 - Buffer volumes
- Remote locations
- Limited space/accessibility
- Climatic conditions sub-optimal
- Standalone operation > 1 month
- Multiple gases



Winderlich et al, 2010. doi:10.5194/amt-3-1113-2010, 2010

Instrument development: Ecotech Spectronus in-situ FTIR

- Upper section**
Computer and temperature control
- Perspex Shipping covers**
(must be removed before operation)
- Middle section**
Optical system – temperature controlled.
- Lower section**
Gas handling system





120 kg
 pump 17 kg
 1160x486x885 mm
 200 W
 Pump 310 W
 24 m multipass cell
 3.5 liter **glass** cell
 0.5-1.5 l.min⁻¹ or static
 MCT detector
 2000-7800 cm⁻¹
 Spectral res 1.0 cm⁻¹
 15 s-60 min averaging
 10-40 °C oper. temp
 Built-in sample drying


Gases	Precision (5 min average)
CO ₂ /ppmv	0.04
CH ₄ /ppbv	0.2
N ₂ O/ppbv	< 0.06
CO/ppbv	0.2
δ ¹³ C-CO ₂ /‰	< 0.08*
δD in H ₂ O vapour/‰	< 1

Improvements tested at ECN


- Aluminium Cell (Bruker)
 - Smaller cell volume 2.5 liter
 - Improved thermal properties
- Cell temperature control
 - Heating at 6 locations, 3 individually controlled zones
 - Increase accuracy of temp meas. with thermistors (0.1- \rightarrow 0.001 K)
 - Temperature control of cell within 10 mK in 60 sec
 - Enhanced air circulation in measurement compartment
- Polishing for reduced air-cell active surface area
- Surface gold plating to reduce wall interactions
- Additional material costs less than € 1000



ICOS |  **INTEGRATED
CARBON
OBSERVATION
SYSTEM**




2015 InGOS Science Conference








2015 InGOS

ICOS |  INTEGRATED
CARBON
OBSERVATION
SYSTEM





Spectronus static mode Allan variance analysis 24 Jan
N2O ppb








ECN, ATV 2014/01/24





Spectronus static mode Allan variance analysis 23 Apr
N2O ppb








ECN, ATV 2014/04/26





Spectronus flow mode Allan variance analysis 24 Apr
N2O ppb







ECN, ATV 2014/04/26



Spectronus metal cell, increased temp control

Species	Allan Var				Precision stdev 3 days	Drift per day	Unit
	flow 1 min	flow 5 min	static 1 min	static 5 min			
CO ₂	0.018	0.007	0.018	0.014	0.031	0.021	ppm
13CO ₂	0.08	0.04	0.03	0.02	0.07	0.03	permille
CH ₄	0.18	0.10	0.20	0.10	0.18	0.11	ppb
N ₂ O	0.15	0.07	0.12	0.05	0.08	0.009	ppb
CO	0.25	0.12	0.20	0.07	0.14	0.04	ppb



Buffer volumes

Two controlling programs


- Temperature control program, three zones, 10 mK accuracy
 - Controls ICP DAS 7005 unit
 - Reads thermistors
 - Switches heating pads (3 W)
 - Feedbacks average cell temperature to Analog Input of Spectronus
- Cycle steering program
 - Switches Valco 12 multiport valve at serial port
 - Triggers Spectronus cycle by pulse of $2V + \text{port number} * 0.2 V$
 - Retrieves molar ratios from Spectronus database
 - Performs at NRT target correction and smoothed calibration
 - Planned: reads Picarro data and sets valve id



- 5 minute cycle
 - Evacuate
 - Flush cell to 500 mbar
 - Evacuate
 - Fill cell to 1100 mbar
 - Collect spectra for 1.5 minute
- 3 full gradients at 4 heights per hour! (4 or 6 at 3 resp. 2 heights)
- Gas use per cycle: 4 liter
- 2000 calibrations/target measurements per 50l tank
- Lifetimes tanks
 - Hi freq target (12h): 2-3 years
 - Lo freq target (1/week): 20-40 years
 - Calibration (2/week): 10 years
- Virtual air sample archive with full IR spectra







Buffer volume
19.5 l
Stainless steel
€ 170,= piece
Max 6 bar










Picarro G2301 vs Spectronus, CBW 200 m




CO2 Cabauw 20m Picarro vs. Spectronus




PROVISIONAL DATA

Cabauw CO2 200 m comparison Spectronus vs. Picarro



Cabauw CH4 200 m comparison Spectronus vs. Picarro



PROVISIONAL DATA



Summary



- Performance improved considerably
- WMO compatibility targets achieved
- Robust deployment at tall tower
- Low maintenance (once per month)
- Very low gas consumption
- Manufacturer will adopt changes (upgrade possible)
- Virtual air archive!