

Real-time analysis of $\delta^{13}\text{C}$ - and $\delta\text{D-CH}_4$ in ambient air with laser spectroscopy: Method development and first Intercomparison results



Simon Eyer, Béla Tuzson, Eliza Harris, Lukas Emmenegger and Joachim Mohn, Empa
Elena Popa, Carina van der Veen, Thomas Röckmann, Utrecht University
Rebecca Fisher, David Lowry, Euan G. Nisbet, RHUL
Willi A. Brand, MPI for Biogeochemistry
Matthias S. Brennwald, Eawag
Hubertus Fischer, University of Bern




Universiteit Utrecht




eawag
aquatic research **ooo**


The global CH₄ cycle



CH₄ Isotopes






InGOS Goals:
0.1 ‰ for $\delta^{13}\text{C}-\text{CH}_4$
1 ‰ for $\delta\text{D}-\text{CH}_4$




Precision \approx ppt (10^{-12})


Analytical Technique





CH₄-Preconcentration





CH₄-Separation





Laser spectroscopy



- two cw-DFB-QCLs @ 7.7 μm
- MPC of 76 m optical path
- TEC MCT IR-detector



Laser spectroscopy – Allan precision



Compact & automated measurement system



QCLAS

Electronics (QCLAS)

TRES

Coolers/Pumps

- CH_4 , $\delta^{13}\text{C}$ & δD
- Compact
- Automated (LabView)
- ~1 hour time resolution

Field Campaign in Dübendorf



Field Campaign in Dübendorf



- TREX-QCLAS (Empa)
- Sampling + IRMS (UU, RHUL, MPI)
- CRDS, OA-ICOS




Universiteit Utrecht


eawag
aquatic research ooo





Repeatability



Temporal trends




Compatibility



--- Extended WMO compatibility goals

Source processes



Conclusion & Outlook



- ✓ Successful CH₄ preconcentration by TREX
- ✓ High precision CH₄ isotopic analysis by QCLAS
≈0.1 ‰ for δ¹³C and <0.5 ‰ for δD-CH₄
- ✓ Repeatability of TREX + QCLAS
0.19 ‰ for δ¹³C and 1.9 ‰ for δD-CH₄
- ✓ Good compatibility of TREX+QCLAS with IRMS for δ¹³C-CH₄ and δD
- ✓ Successful determination of source processes
- ✓ 6 months field campaign at Cabauw, NL

Special thanks to...

- Empa
 - J. Mohn, B. Tuzson und L. Emmenegger
 - S. Reimann, M. Vollmer, M. Hill
 - C. Zellweger and M. Steinbacher, GAW
 - D. Brunner, Eliza Harris
 - Hydrogen Group, Workshop
- University of Utrecht
 - Thomas Röckmann, Carina van der Veen, Elena Popa
- MPI Jena
 - Willi A. Brand
- RHUL
 - Rebecca Fisher, David Lowry, Euan G. Nisbet
- University Bern
 - H. Fischer
- INGOS & SNF



Universiteit Utrecht



eawag
aquatic research ooo

Special thanks to...



- Empa
 - J. Mohn, B. Tuzson und L. Emmenegger
 - S. Reimann, M. Vollmer, M. Hill
 - C. Zellweger and M. Steinbacher, GAW
 - D. Brunner, Eliza Harris
 - Hydrogen Group, Workshop
- University of Utrecht
 - Thomas Röckmann, Carina van der Veen, Elena Popa
- MPI Jena
 - Willi A. Brand
- RHUL
 - Rebecca Fisher, David Lowry, Euan G. Nisbet
- University Bern
 - H. Fischer
- INGOS & SNF




Thanks for listening!!






Universiteit Utrecht




Setup



Compatibility – preliminary results



Compact & automated measurement system




Participants campaign Dübendorf


- $\delta^{13}\text{C}$ and δD and CH_4 by QCLAS
 CH_4 by CRDS (Picarro)
- $\delta^{13}\text{C}$ and CH_4 by OA-ICOS (LGR)
 $\delta^{13}\text{C}$ and δD and CH_4 by flask sampling + IRMS
- $\delta^{13}\text{C}$ and δD and CH_4 by flask sampling + IRMS
- $\delta^{13}\text{C}$ and CH_4 by bag sampling + IRMS
- $\delta^{13}\text{C}$ and CH_4 by CRDS (Picarro)




The global CH₄ cycle




CH₄ Isotopes




98.9%



1.1%



0.06%



Cabauw - Messkampagne



- Oct. 2014 – March 2015
- 20 m height

Compact & automated measurement system



Optics

Electronics (QCLAS)

TREX

Coolers/Pumps


- $\delta^{13}\text{C}$ & δD
- Kompakt
- Vollautomatisch (LabView)
- <2.4kW
- ≈ 350 kChf

Summary and Outlook




- ✓ Successful CH₄ preconcentration by TREX
- ✓ High precision CH₄ isotopic analysis by QCLAS
≈0.1 ‰ for δ¹³C and <0.5 ‰ for δD-CH₄
- ✓ Repeatability of TREX + QCLAS
0.19 ‰ for δ¹³C and 1.9 ‰ for δD-CH₄
- ✓ 4 months field campaign at Cabauw, NL
- ✓ Good compatibility of TREX+QCLAS with IRMS for δ¹³C-CH₄ and δD
- ✓ Mainly ruminant sources

Temporal trends




Compatibility – preliminary results



Overview

- Introduction
- Analytical Technique
- Campaign in Dübendorf
- Conclusion and Outlook



CH₄-Extraktion

